



Ready-to-use wipe for alcohol-free surface disinfection and cleaning

- Residual antimicrobial effect for up to 10 days
- effective against Coronavirus
- excellent skin tolerance
- high cleaning effect
- practical applicability tested in the 4-field-test



Intended use

Bacoban^{DL} wipe is a ready-to-use, waterbased wipe for surface disinfection and cleaning with residual, antimicrobial effect.

Directions for use

Use the **Bacoban^{DL}** wipe on the surface, rub / spread evenly, then allow to dry. Do not rinse or wipe off. In case of heavy soiling, use several wipes

Residual antimicrobial effect for up to 10 days (tested according to ASTM E 2180).
Dermatologically tested : excellent skin tolerance.

Ingredients

100 g of solution containing: Benzalkonium chloride 0,26 g, sodium pyrithione 0,025 g, polycondensates, perfume substances, purified water.

Ingredients according to regulation (EC) No 648/2004: less than 5 % cationic surfactants, disinfectants, perfumes. Does not contain formaldehyde. Does not contain phenols.

Storage

Store **Bacoban^{DL}** wipe at +5 bis +25°C. Protect from frost and direct sunlight.

Test reports

Our products are tested extensively. We are pleased to provide you with our microbiological and dermatological reports on request.

Safety notice

Safety data sheet available on request.

Important! The residual antimicrobial effect of **Bacoban^{DL}** is no excuse for not carrying out necessary disinfection measures in cases of visible and/or questionable contamination!

Use disinfectants safely. Always read label and product information before use.

DISTRIBUTOR

Efficacy

Testing according to EN14476 and guideline DVV / RKI	time of exposure
Virucidal activity against enveloped viruses (incl. Coronavirus , HIV, HBV, HCV, Influenza, Herpes, Ebolavirus...)	5 min
Rotavirus	5 min
Norovirus	30 min
Human papillomavirus (HPV)	60 min
VAH list (tested according to VAH-guidelines 2015)	
bactericidal and yeasticidal	5 min
EN guidelines	
EN 1040 (bactericidal)	5 min
EN 1275 (yeasticidal)	5 min
EN 1276 (bactericidal)	5 min
EN 13727 (bactericidal)	5 min
EN 13624 (yeasticidal, fungicidal)	5 min
EN 1650 (yeasticidal)	5 min
EN 13697 (bactericidal)	15 min
EN 14561 (bactericidal)	60 min
EN 14562 (yeasticidal)	60 min
EN 16615 (bactericidal, yeasticidal) (4 field test)	5 min